



225-273-1900 • Fax: 225-273-5555
 info@stillmeinc.com

Patient: Last Name: _____ First Name _____

Fitter: _____ Clinic _____


Fitter Title: _____ (example: PT/OT/PTA)

Date: _____

JOBST® FarrowWrap® Basic

3 Easy Measuring Steps:

1. Length (A-D) - Floor to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) - Calf (widest)
3. Circumference (B) - Ankle (smallest part)



Measurements		
	Left	Right
cC		
cB		
ℓA-D		

JOBST® FarrowWrap® Basic (30-40 mmHg*)

JOBST® FarrowWrap® Basic Size Chart

	XS	S	M	L	XL
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
A-D Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm

Ref No.	Size	Length	Quantity
7666000	XS	Regular	
7666001	S	Regular	
7666002	M	Regular	
7666003	L	Regular	
7666004	XL	Regular	
7666005	XS	Tall	
7666006	S	Tall	
7666007	M	Tall	